SPANISH

Para utilizar con el TearLab Osmolarity System	IVD	REF Each level sold separately
	Uso diagnóstico in vitro	Normal - 100012 Alto - 100013

USO PREVISTO

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Authorized Representative for the
European Union:
Cavendish Scott Ltd
PO Box 107
SG5 1FW, England

One or more of the following patents
may apply: U.S. Patents 7,017,394; 7,051,569;
7,111,502; 7,129,717; 7,204,122

ProClin® is a registered trademark of
Rohm and Haas Company.

Authorized Representative

Catalog Number

Consult Instructions for Use

Control

European Community

For In Vitro Diagnostic Use

Manufacturer

Store At

ENGLISH

For use with the TearLab Osmolarity System	IVD	REF Each level sold separately
	For In Vitro Diagnostic Use	Normal - 100012 High - 100013

INTENDED USE

TearLab Osmolarity Control Solution is a quality control material used for verifying the quantitative performance of TearLab Osmolarity Test Cards when used with the TearLab Osmolarity System.

Routine use of TearLab Osmolarity Control Solution enables laboratories to monitor day-to-day test variation and lot-to-lot test kit performance and to assist in operator training. These controls can also be used to troubleshoot invalid results and to identify increases in random or systemic error.

PRODUCT DESCRIPTION

TearLab Osmolarity Control Solution is provided in two (2) levels, either a Normal Level kit (Product Number 100012) or a High Level kit (Product Number 100013), for monitoring Test Card performance at the normal and high clinical range of tear osmolarity. It is packaged in sealed glass ampoules, each containing 1.0 mL of solution. These ampoules are intended for onetime use. Discard any solution that remains after use.

ACTIVE INGREDIENTS

This product is an aqueous solution containing salt (NaCl) traceable to National Institute of Standards and Technology (NIST) and ProClin® 300 preservative and contains no human or biological materials.

STORAGE

↕^{30°} The expiration date on the package is for product stored at a room temperature of 2°C-30°C (36°F-86°F). Avoid freezing.

DIRECTIONS FOR USE

- Allow the control solutions to come to room temperature.
- Tap the neck of the control solution ampoule to remove any fluid. Using the blue rubber “sleeve,” snap off the neck of the ampoule and discard in a sharps container. Retain the sleeve for future use.
- Attach a TearLab Osmolarity Test Card onto a TearLab Pen. A green light on top of the Pen will illuminate and the Pen will beep.
- Invert the ampoule – the control fluid will not spill out. Remove the protective cover from the Test Card and dip the tip of the Pen into the control fluid until it beeps and the green light turns off.
- Dock the Pen into the TearLab Reader. Confirm that the test result is within the expected range indicated below.

EXPECTED VALUES

Listed below are the nominal values and the corresponding expected ranges for the control materials, assuming nominal performance by the TearLab Osmolarity System.

	EXPECTED (TARGET)	VALUES EXPECTED RANGES
Normal	297 mOsms/L	±15 mOsms/L (±3SD)
High	338 mOsms/L	±15 mOsms/L (±3SD)

LIMITATIONS

- This product is intended for use as a quality control material in evaluating the performance of the TearLab Osmolarity System. It is not suitable for use on other osmometers that do not use TearLab Osmolarity Test Cards for osmolarity determinations.
- TearLab Osmolarity Control Solutions should complement other aspects of a complete quality control program.
- TearLab Osmolarity Control Solutions are intended for use in undiluted form. Do not dilute the solutions prior to use.
- TearLab Osmolarity Control Solutions are intended for single use only. Discard solution after use.
- Use of expired control solutions may produce erroneous results. Do not use beyond the expiration date on the package.

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©2009 TearLab Corporation