For ordering information, contact your local sales representative or TearLab Customer Support in your area.

#### OSMOLARITY REFERENCES

1. Pederseen-Bjergaard K, Smidt BC, "Electrical conductivity, osmotic pressure, and hydrogen ion concentration of human lachrymal fluid," Acta Derm Venereol Suppl (Stockh) 1952; 32(29): 261-7.
2. American Academy of Ophthalmology, Information Statement, "Minimizing Transmission of Blood-borne Pathogens and Surface Infectious Agents in Ophthalmic Offices and Operating Rooms," March 1, 2002 (www.aaoo.org).
3. Tomlinson A. Tear Film Osmolarity: Determination of a Referent for Dry Eye Diagnosis. IOVS 47(10), 2006.

 TearLab Corporation, San Diego, CA, USA, 1-858-455-6006, [www.tearlab.com](http://www.tearlab.com)

 Cavendish Scott LTD, PO Box 107, SG5 1FW, England

 For complete instructions for use, please refer to the TearLab User Manual

CE

US Pat. # 7017394, 7051569, 7111502, 7129717, 7204122

 TearLab™

930018REV C

©2009 TearLab Corporation

TearLab™

## TEST CARDS ENGLISH

*For professional In-Vitro diagnostic use.*

### INTENDED USE

The TearLab™ Osmolarity Test is a tear fluid collection and testing device utilizing nanofluidic lab-on-a-chip technology for the measurement of osmolarity (concentration of dissolved, active particles in solution) of human tears in normal and dry eye patients.

### PRINCIPLES OF THE PROCEDURE

The TearLab Osmolarity test utilizes an impedance measurement of tear fluid to provide a calculated measurement of osmolarity.<sup>1</sup>

### MATERIALS

TearLab™ Osmolarity Test Card: a single-use, individually packaged, non-sterile, microchip to enable measurement of the impedance of the tear fluid sample. Each test card is hygienically clean and contains a protective cover. Test Cards are designed to work in conjunction with the TearLab Osmolarity System.

### REAGENTS

Test Cards contain no reagents or chemicals.

### MATERIALS NOT PROVIDED

TearLab Osmolarity System consisting of: Reader, Pens, and Electronic Test Cards Normal and High Osmolarity Quality Control Solutions

### PRECAUTIONS

- Test Cards are single use only.
- Use only at ambient temperature of 15 – 30 °C (59 – 86 °F).
- Proper handling and disposal methods of used Test Cards should be established according to relevant state and federal regulations.
- TearLab is designed to collect samples of tear fluid from the eye, a non-sterile environment. The AAO has issued guidance to minimize the transmission of ocular surface infectious agents.<sup>2</sup> Prevention of transmission of these pathogens requires good hygienic techniques, such as hand washing and the cleanliness of instruments that would contact the eye. TearLab Osmolarity Test Cards are single-use, hygienically clean, contain a protective cover and are individually packaged. No cleaning is required.

### STORAGE AND STABILITY

Store Osmolarity Test Cards at room temperature 15-30°C (59-86°F).

**IMPORTANT:** Test Cards not stored at room temperature should be at ambient temperature before testing. Cold Test Cards may require several minutes to reach ambient temperature. Osmolarity Test Cards are stable until the expiration date marked on the label.

### SAMPLE MATERIAL

Human tear fluid samples may be used. Collect tear fluid samples directly from the eye. Do not collect tear fluid samples from patients that have used eye drops within 2 hours prior to testing.

Do not collect or store tear fluid samples for transport or testing at a later time.

Do not collect tear fluid after ocular surface staining.

Do not collect tear fluid after other invasive ocular diagnostic testing.

Do not collect tear fluid within ten minutes after a slit lamp examination.

Do not collect tear fluid from a patient that has been crying.

**TESTING PROCEDURE** (Refer to the TearLab Osmolarity System User Manual for more information)

### QUALITY CONTROL

Electronic Check Card: Refer to the TearLab Osmolarity System User Manual for more information.

**OSMOLARITY CONTROL SOLUTIONS:** Refer to TearLab Osmolarity System User Manual for more information.

If correct electronic or control solution test results are not obtained, do not test patient samples and contact your local sales representative or TearLab Customer Support in your area.

### EXPECTED RESULTS

TearLab measurement range is linear from 275 mOsm/L to 400 mOsm/L. Test results outside of this measurement range will be reported as either “Below Range” indicating a measurement below 275 mOsm/L or “Above Range” indicating a measurement above 400 mOsm/L.

Reference tear osmolarity values for normal and dry eye disease patients:

- Normal 275 mOsm/L to 316 mOsm/L (mean 302 mOsm/L)<sup>3</sup>
- Dry Eye Disease > 316 mOsm/L (mean 327 mOsm/L)<sup>3</sup>

### PRECISION *Within-run, Between-day and Total Precision*

Precision was determined over 10 days with three different osmolarity levels of contrived tear samples. Two times two replicates were analyzed per day of each level in accordance with the Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2.

	WITHIN RUN	DAY-TO-DAY	TOTAL
<b>296 mOsm/L</b>			
SD (mOsm/L)	3.7	3.4	7.1
CV (%)	1.26%	1.14%	2.40%
<b>315 mOsm/L</b>			
SD (mOsm/L)	4.6	1.9	6.5
CV (%)	1.47%	0.61%	2.07%
<b>360 mOsm/L</b>			
SD (mOsm/L)	4.4	1.1	5.5
CV (%)	1.21%	0.32%	1.53%

### ACCURACY (Method Comparison)

The correlation study was performed using 28 contrived tear samples of various osmolarity levels within the clinical reference range on both the TearLab Osmolarity System and the Wescor Model 5520 vapor pressure osmometer calibrated to National Institute of Standards Institute (NIST) traceable standards.

#### RESULTS

Range of values (mOsm/L)	Number of Samples	Slope	Y Intercept	Correlation Coefficient
<b>280-370</b>	<b>28</b>	<b>0.95</b>	<b>10.69</b>	<b>0.98</b>

#### INTERFERING SUBSTANCES

Tear samples are to be collected directly from the eye. Do not collect tears within two hours of use of eye drops. Do not collect tears after ocular surface staining. Do not collect tears from patients who have been crying.

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